

BDKRB1 Antibody



PACO04973

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

WB:1:500-1:2000, IF:1:200-1:1000

Protein Background:**Gene ID:**

BDKRB1

Uniprot

P46663

Synonyms:

BDKRB1; BRADYB1; B1 bradykinin receptor; B1R; BK-1 receptor

Immunogen:

Synthesized peptide derived from the Internal region of human Bradykinin B1 R.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A