

SAA4 Antibody



PACO04984

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

IHC:1:100-1:300, IHC:1:200-1:1000

Protein Background:**Gene ID:**

SAA4

Uniprot

P35542

Synonyms:

SAA4; CSAA; Serum amyloid A-4 protein; Constitutively expressed serum amyloid A protein; C-SAA

Immunogen:

Synthesized peptide derived from the C-terminal region of human SAA4.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A