

Phospho-ATXN1 (S776) Antibody



PACO04985

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300,
IF:1:200-1:1000

Protein Background:

Gene ID:

ATXN1

Uniprot

P54253

Synonyms:

ATXN1; ATX1; SCA1; Ataxin-1; Spinocerebellar ataxia type 1 protein

Immunogen:

Synthesized peptide derived from human Ataxin-1 around the phosphorylation site of S776.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A