

Phospho-BID (S78) Antibody



PACO04997

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

IHC:1:100-1:300

Protein Background:

Gene ID:

BID

Uniprot

P55957

Synonyms:

BID; BH3-interacting domain death agonist; p22 BID; BID

Immunogen:

Synthesized peptide derived from human BID around the phosphorylation site of S78.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A