

SRR Antibody



PACO05000

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:5000, IHC:1:100-1:300

Protein Background:**Gene ID:**

SRR

Uniprot

Q9GZT4

Synonyms:

SRR; Serine racemase; D-serine ammonia-lyase; D-serine dehydratase; L-serine ammonia-lyase; L-serine dehydratase

Immunogen:

Synthesized peptide derived from the Internal region of human Serine racemase.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A