

Phospho-SDC4 (S179) Antibody



PACO05002

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:20000, IHC:1:100-1:300, IF:1:200-1:1000

Protein Background:

Gene ID:

SDC4

Uniprot

P31431

Synonyms:

SDC4; Syndecan-4; SYND4; Amphiglycan; Ryudocan core protein

Immunogen:

Synthesized peptide derived from human Syndecan-4 around the phosphorylation site of S179.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A