

ABHD4 Antibody



PACO05003

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:20000, IHC:1:100-1:300

Protein Background:**Gene ID:**

ABHD4

Uniprot

Q8TB40

Synonyms:

ABHD4; Abhydrolase domain-containing protein 4; alpha/beta-hydrolase 4; Lyso-N-acylphosphatidylethanolamine lipase

Immunogen:

Synthesized peptide derived from the C-terminal region of human ABHD4.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A