

# SRSF5 Antibody



PACO05023

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC, IF

**Recommended dilutions:**

ELISA:1:20000, IHC:1:100-1:300, IF:1:200-1:1000

**Protein Background:****Gene ID:**

SRSF5

**Uniprot**

Q13243

**Synonyms:**

SRSF5; HRS; SFRS5; SRP40; Serine/arginine-rich splicing factor 5; Delayed-early protein HRS; Pre-mRNA-splicing factor SRP40; Splicing factor; arginine/serine-rich 5

**Immunogen:**

Synthesized peptide derived from the Internal region of human SRp40.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A