

ARHGAP9 Antibody



PACO05028

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:40000, WB:1:500-1:2000,
IHC:1:100-1:300

Protein Background:

Gene ID:

ARHGAP9

Uniprot

Q9BRR9

Synonyms:

ARHGAP9; Rho GTPase-activating protein 9; Rho-type GTPase-activating protein 9

Immunogen:

Synthesized peptide derived from the Internal region of human ARHGAP9.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A