

NOM1 Antibody



PACO05039

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:40000, IHC:1:100-1:300

Protein Background:

Gene ID:

NOM1

Uniprot

Q5C9Z4

Synonyms:

NOM1; C7orf3; Nucleolar MIF4G domain-containing protein 1; SGD1 homolog

Immunogen:

Synthesized peptide derived from the C-terminal region of human NOM1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A