

# GCNT6 Antibody



PACO05040

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IF

**Recommended dilutions:**

ELISA:1:40000, WB:1:500-1:2000, IF:1:200-1:1000

**Protein Background:**

**Gene ID:**

GCNT6

**Uniprot**

Q5T4J0

**Synonyms:**

GCNT6; beta-1; 3-galactosyl-O-glycosyl-glycoprotein beta-1,6-N-acetylglucosaminyltransferase 6

**Immunogen:**

Synthesized peptide derived from the Internal region of human GCNT6.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A