

# NKX2-4 Antibody



PACO05042

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Monkey

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

ELISA:1:40000, WB:1:500-1:2000

**Protein Background:**

**Gene ID:**

NKX2-4

**Uniprot**

Q9H2Z4

**Synonyms:**

NKX2-4; NKX2D; Homeobox protein Nkx-2.4; Homeobox protein NK-2 homolog D

**Immunogen:**

Synthesized peptide derived from the C-terminal region of human NKX2.4.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A