

WAVE4 Antibody



PACO05044

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:10000, IHC:1:100-1:300

Protein Background:

Gene ID:

WAVE4

Uniprot

Q8IV90

Synonyms:

Putative Wiskott-Aldrich syndrome protein family member 4

Immunogen:

Synthesized peptide derived from the Internal region of human WAVE4.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A