

SHB Antibody



PACO05048

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:5000, WB:1:500-1:2000

Protein Background:**Gene ID:**

SHB

Uniprot

Q15464

Synonyms:

SHB; SH2 domain-containing adapter protein B

Immunogen:

Synthesized peptide derived from the C-terminal region of human Shb.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A