

# Phospho-SHC1 (S36) Antibody



PACO05050

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:40000, IHC:1:100-1:300

**Protein Background:****Gene ID:**

SHC1

**Uniprot**

P29353

**Synonyms:**

SHC1; SHC; SHCA; SHC-transforming protein 1; SHC-transforming protein 3; SHC-transforming protein A; Src homology 2 domain-containing-transforming protein C1; SH2 domain protein C1

**Immunogen:**

Synthesized peptide derived from human Shc around the phosphorylation site of S36.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A