

# Phospho-ANAPC1 (S355) Antibody



PACO05052

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:5000, IHC:1:100-1:300

**Protein Background:****Gene ID:**

ANAPC1

**Uniprot**

Q9H1A4

**Synonyms:**

ANAPC1; TSG24; Anaphase-promoting complex subunit 1; APC1; Cyclosome subunit 1; Mitotic checkpoint regulator; Testis-specific gene 24 protein

**Immunogen:**

Synthesized peptide derived from human APC1 around the phosphorylation site of S355.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A