

# IFT22 Antibody



PACO05059

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:5000, IHC:1:100-1:300

**Protein Background:**

**Gene ID:**

IFT22

**Uniprot**

Q9H7X7

**Synonyms:**

RABL5; Rab-like protein 5

**Immunogen:**

Synthesized peptide derived from the C-terminal region of human Rab L5.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A