

SIX1 Antibody



PACO05066

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:20000, WB:1:500-1:2000

Protein Background:**Gene ID:**

SIX1

Uniprot

Q15475

Synonyms:

SIX1; Homeobox protein SIX1; Sine oculis homeobox homolog 1

Immunogen:

Synthesized peptide derived from the Internal region of human Six1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A