

MRPS24 Antibody



PACO05067

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:40000, IHC:1:100-1:300

Protein Background:**Gene ID:**

MRPS24

Uniprot

Q96EL2

Synonyms:

MRPS24; HSPC335; 28S ribosomal protein S24; mitochondrial; MRP-S24; S24mt; bMRP-47; bMRP47

Immunogen:

Synthesized peptide derived from the Internal region of human MRP-S24.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A