

MRPS5 Antibody



PACO05068

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:40000, IHC:1:100-1:300

Protein Background:**Gene ID:**

MRPS5

Uniprot

P82675

Synonyms:

MRPS5; 28S ribosomal protein S5; mitochondrial; MRP-S5; S5mt

Immunogen:

Synthesized peptide derived from the C-terminal region of human MRP-S5.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A