

SLC6A6 Antibody



PACO05077

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:40000, WB:1:500-1:2000,
IHC:1:100-1:300

Protein Background:**Gene ID:**

SLC6A6

Uniprot

P31641

Synonyms:

SLC6A6; Sodium- and chloride-dependent taurine transporter; Solute carrier family 6 member 6

Immunogen:

Synthesized peptide derived from the C-terminal region of human TAUT.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A