

SLC6A11 Antibody



PACO05084

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:10000, WB:1:500-1:2000

Protein Background:**Gene ID:**

SLC6A11

Uniprot

P48066

Synonyms:

SLC6A11; GABT3; GAT3; Sodium- and chloride-dependent GABA transporter 3; GAT-3; Solute carrier family 6 member 11

Immunogen:

Synthesized peptide derived from the Internal region of human GABA T-3.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A