

# BMP6 Antibody



PACO05085

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

ELISA:1:40000, WB:1:500-1:2000,  
IHC:1:100-1:300

**Protein Background:****Gene ID:**

BMP6

**Uniprot**

P22004

**Synonyms:**

BMP6; VGR; Bone morphogenetic protein 6; BMP-6; VG-1-related protein; VG-1-R;  
VGR-1

**Immunogen:**

Synthesized peptide derived from the Internal region of human BMP-6.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A