

SOS2 Antibody



PACO05105

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:10000, WB:1:500-1:2000

Protein Background:**Gene ID:**

SOS2

Uniprot

Q07890

Synonyms:

SOS2; Son of sevenless homolog 2; SOS-2

Immunogen:

Synthesized peptide derived from the Internal region of human Sos 2.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A