

Phospho-SOX9 (S181) Antibody



PACO05106

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

ELISA:1:5000, WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000

Protein Background:

Gene ID:

SOX9

Uniprot

P48436

Synonyms:

SOX9; Transcription factor SOX-9

Immunogen:

Synthesized peptide derived from human Sox-9 around the phosphorylation site of S181.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A