

SOX12 Antibody



PACO05107

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:20000, IHC:1:100-1:300

Protein Background:**Gene ID:**

SOX12

Uniprot

O15370

Synonyms:

SOX12; SOX22; Transcription factor SOX-12; Protein SOX-22

Immunogen:

Synthesized peptide derived from the Internal region of human Sox-12.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A