

Phospho-FOXL2 (S263) Antibody



PACO05112

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:20000, WB:1:500-1:2000

Protein Background:

Gene ID:

FOXL2

Uniprot

P58012

Synonyms:

FOXL2; Forkhead box protein L2

Immunogen:

Synthesized peptide derived from human FoxL2 around the phosphorylation site of S263.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A