

Phospho-BRCA1 (S1423) Antibody



PACO05115

Product Information

Size:

50ug

Reactivity:

Human, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:5000, WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:

Gene ID:

BRCA1

Uniprot

P38398

Synonyms:

BRCA1; RNF53; Breast cancer type 1 susceptibility protein; RING finger protein 53

Immunogen:

Synthesized peptide derived from human BRCA1 around the phosphorylation site of S1423.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A