

BRCA1 Antibody



PACO05118

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:20000, IHC:1:100-1:300

Protein Background:

Gene ID:

BRCA1

Uniprot

P38398

Synonyms:

BRCA1; RNF53; Breast cancer type 1 susceptibility protein; RING finger protein 53

Immunogen:

Synthesized peptide derived from human BRCA1 around the non-phosphorylation site of S988.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A