

Phospho-SRF (S103) Antibody



PACO05120

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:5000, WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:**Gene ID:**

SRF

Uniprot

P11831

Synonyms:

SRF; Serum response factor; SRF

Immunogen:

Synthesized peptide derived from human SRF around the phosphorylation site of S103.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A