

# Phospho-BRAF (T599) Antibody



PACO05124

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:40000, IHC:1:100-1:300

**Protein Background:****Gene ID:**

BRAF

**Uniprot**

P15056

**Synonyms:**

BRAF; BRAF1; RAFB1; Serine/threonine-protein kinase B-raf; Proto-oncogene B-Raf; p94; v-Raf murine sarcoma viral oncogene homolog B1

**Immunogen:**

Synthesized peptide derived from human Raf-B around the phosphorylation site of T599.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A