

PACO05128

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:20000, WB:1:500-1:2000,
IHC:1:100-1:300

Protein Background:

Gene ID:

SRY

Uniprot

Q05066

Synonyms:

SRY; TDF; Sex-determining region Y protein; Testis-determining factor

Immunogen:

Synthesized peptide derived from the Internal region of human SRY.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A