

# Phospho-STAT5A/STAT5B (Y694/699) Antibody



PACO05141

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

ELISA:1:20000, WB:1:500-1:2000,  
IHC:1:100-1:300

**Protein Background:****Gene ID:**

STAT5A/STAT5B

**Uniprot**

P42229/P51692

**Synonyms:**

STAT5A; STAT5; Signal transducer and activator of transcription 5A; STAT5B; Signal transducer and activator of transcription 5B

**Immunogen:**

Synthesized peptide derived from human Stat5 around the phosphorylation site of Y694/699.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A