

BRS3 Antibody



PACO05149

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

ELISA:1:10000, WB:1:500-1:2000,
IHC:1:100-1:300, IF:1:200-1:1000

Protein Background:**Gene ID:**

BRS3

Uniprot

P32247

Synonyms:

BRS3; Bombesin receptor subtype-3; BRS-3

Immunogen:

Synthesized peptide derived from the Internal region of human BRS-3.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A