

Phospho-SYK (Y323) Antibody



PACO05154

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

ELISA:1:5000, WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000

Protein Background:**Gene ID:**

SYK

Uniprot

P43405

Synonyms:

SYK; Tyrosine-protein kinase SYK; Spleen tyrosine kinase; p72-Syk

Immunogen:

Synthesized peptide derived from human Syk around the phosphorylation site of Y323.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A