

# Phospho-SYN1 (S62) Antibody



PACO05158

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

ELISA:1:10000, WB:1:500-1:2000,  
IHC:1:100-1:300

**Protein Background:**

**Gene ID:**

SYN1

**Uniprot**

P17600

**Synonyms:**

SYN1; Synapsin-1; Brain protein 4.1; Synapsin I

**Immunogen:**

Synthesized peptide derived from human Synapsin I around the phosphorylation site of S62.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A