

# TAZ Antibody



PACO05175

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

ELISA:1:5000, WB:1:500-1:2000

**Protein Background:**

**Gene ID:**

TAZ

**Uniprot**

Q16635

**Synonyms:**

TAZ; EFE2; G4.5; Tafazzin; Protein G4.5

**Immunogen:**

Synthesized peptide derived from the Internal region of human Tafazzin.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A