

TESK1 Antibody



PACO05198

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:40000, WB:1:500-1:2000

Protein Background:**Gene ID:**

TESK1

Uniprot

Q15569

Synonyms:

TESK1; Dual specificity testis-specific protein kinase 1; Testicular protein kinase 1

Immunogen:

Synthesized peptide derived from the Internal region of human TESK1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A