

# Phospho-TGFBR3 (T842) Antibody



PACO05201

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC, IF

**Recommended dilutions:**

ELISA:1:5000, IHC:1:100-1:300, IF:1:200-1:1000

**Protein Background:****Gene ID:**

TGFBR3

**Uniprot**

Q03167

**Synonyms:**

TGFBR3; Transforming growth factor beta receptor type 3; TGF-beta receptor type 3; TGFR-3; betaglycan; Transforming growth factor beta receptor III; TGF-beta receptor type III

**Immunogen:**

Synthesized peptide derived from human TGF & beta; RIII around the phosphorylation site of T842.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A