

TGFBR3 Antibody



PACO05202

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IF

Recommended dilutions:

ELISA:1:10000, IF:1:200-1:1000

Protein Background:**Gene ID:**

TGFBR3

Uniprot

Q03167

Synonyms:

TGFBR3; Transforming growth factor beta receptor type 3; TGF-beta receptor type 3; TGF-beta receptor type III; TGF-beta receptor type III

Immunogen:

Synthesized peptide derived from the Internal region of human TGF β RIII.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A