

Phospho-CLDN5 (Y217) Antibody



PACO05216

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:5000, IHC:1:100-1:300

Protein Background:**Gene ID:**

CLDN5

Uniprot

O00501

Synonyms:

CLDN5; AWAL; TMVCF; Claudin-5; Transmembrane protein deleted in VCFS; TMDVCF

Immunogen:

Synthesized peptide derived from human Claudin-5 around the phosphorylation site of Y217.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A