

# TNFAIP2 Antibody



PACO05218

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC, IF

**Recommended dilutions:**

ELISA:1:20000, IHC:1:100-1:300, IF:1:200-1:1000

**Protein Background:**

**Gene ID:**

TNFAIP2

**Uniprot**

Q03169

**Synonyms:**

TNFAIP2; Tumor necrosis factor alpha-induced protein 2; TNF alpha-induced protein 2; Primary response gene B94 protein

**Immunogen:**

Synthesized peptide derived from the Internal region of human TNF &alpha; -IP 2.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A