

# Cleaved-C1R (R463) Antibody



PACO05226

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

ELISA:1:10000, WB:1:500-1:2000

**Protein Background:**

**Gene ID:**

C1R

**Uniprot**

P00736

**Synonyms:**

C1R; Complement C1r subcomponent; Complement component 1 subcomponent r

**Immunogen:**

Synthesized peptide derived from the Internal region of human C1r HC.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A