

Acetyl-TP53 (K386) Antibody



PACO05233

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:40000, IHC:1:100-1:300, IF:1:200-1:1000

Protein Background:**Gene ID:**

TP53

Uniprot

P04637

Synonyms:

TP53; P53; Cellular tumor antigen p53; Antigen NY-CO-13; Phosphoprotein p53; Tumor suppressor p53

Immunogen:

Synthesized peptide derived from human p53 around the acetylation site of K386.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A