

# Phospho-TP53 (T81) Antibody



PACO05236

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:20000, IHC:1:100-1:300

**Protein Background:****Gene ID:**

TP53

**Uniprot**

P04637

**Synonyms:**

TP53; P53; Cellular tumor antigen p53; Antigen NY-CO-13; Phosphoprotein p53; Tumor suppressor p53

**Immunogen:**

Synthesized peptide derived from human p53 around the phosphorylation site of T81.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A