

# Phospho-TP53BP1 (S25) Antibody



PACO05245

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC, IF

**Recommended dilutions:**

ELISA:1:10000, IHC:1:100-1:300, IF:1:200-1:1000

**Protein Background:****Gene ID:**

TP53BP1

**Uniprot**

Q12888

**Synonyms:**

TP53BP1; Tumor suppressor p53-binding protein 1; 53BP1; p53-binding protein 1; p53BP1

**Immunogen:**

Synthesized peptide derived from human 53BP1 around the phosphorylation site of S25.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A