

TRPS1 Antibody



PACO05254

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:5000, IHC:1:100-1:300

Protein Background:**Gene ID:**

TRPS1

Uniprot

Q9UHF7

Synonyms:

TRPS1; Zinc finger transcription factor Trps1; Tricho-rhino-phalangeal syndrome type I protein; Zinc finger protein GC79

Immunogen:

Synthesized peptide derived from the N-terminal region of human TRPS1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A