

# Phospho-TSC2 (Y1571) Antibody



PACO05255

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

ELISA:1:5000, WB:1:500-1:2000, IHC:1:100-1:300

**Protein Background:**

**Gene ID:**

TSC2

**Uniprot**

P49815

**Synonyms:**

TSC2; TSC4; Tuberin; Tuberous sclerosis 2 protein

**Immunogen:**

Synthesized peptide derived from human Tuberin around the phosphorylation site of Y1571.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A