

UMOD Antibody



PACO05281

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:10000, WB:1:500-1:2000

Protein Background:**Gene ID:**

UMOD

Uniprot

P07911

Synonyms:

UMOD; Uromodulin; Tamm-Horsfall urinary glycoprotein; THP

Immunogen:

Synthesized peptide derived from the Internal region of human THP.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A