

# Phospho-VASP (S157) Antibody



PACO05285

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

ELISA:1:5000, WB:1:500-1:2000, IHC:1:100-1:300

**Protein Background:**

**Gene ID:**

VASP

**Uniprot**

P50552

**Synonyms:**

VASP; Vasodilator-stimulated phosphoprotein; VASP

**Immunogen:**

Synthesized peptide derived from human VASP around the phosphorylation site of S157.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A